


ELISA Kits

for Metabolism Diseases Research



//

肥胖症

糖尿病(消渴病)

糖尿病神经病变

糖尿病酮症酸中毒

糖尿病心血管病变

//

Elabscience®

Elabscience Biotechnology Co.,Ltd

代谢疾病

代谢疾病 (代谢病) 是一种影响人类 (或动物) 细胞生产能量的障碍, 是因代谢问题引起的疾病, 包括代谢障碍和代谢旺盛等原因, 又称为新陈代谢失调症。大部分代谢疾病都是遗传性疾病, 而有部分是从饮食、毒素、感染等而有。遗传性的代谢疾病一般称为先天性代谢缺陷。一般来说, 遗传性的代谢疾病都是先天性缺少或不正常构成在细胞代谢过程中重要的酶。


用于代谢疾病研究的 ELISA 试剂盒

Elabscience 根据往期文献将临床 5 种代谢疾病 (肥胖、糖尿病、糖尿病酮症酸中毒等) 研究中的相关指标进行了整理, 便于代谢疾病研究人员能快速查找并获取指标信息。请访问 Elabscience 官网 (www.elabscience.cn), 以获取详细产品信息。

代谢疾病研究中的热点指标:

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-0028c	8-OHdG (8-Hydroxydeoxyguanosine) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Competitive
E-EL-0030c	Cortisol ELISA Kit	12.5--800 ng/mL	7.5 ng/mL	Competitive
E-EL-H0002c	Human ACE(Angiotensin I Converting Enzyme) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0272c	Human ADH(Antidiuretic Hormone) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Competitive
E-EL-H0004c	Human ADP/Acrp30(Adiponectin) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H0295c	Human AGER(Total Advanced Glycosylation End Product Specific Receptor) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0005c	Human AgRP(Agouti Related Protein) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H2181c	Human ALB(Albumin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0329c	Human ANG II R1(Angiotensin II Receptor 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0125c	Human ApoA1(Apolipoprotein A1) ELISA Kit	3.13--200 ng/mL	1.88 ng/mL	Sandwich
E-EL-H0464c	Human ApoB(Apolipoprotein B) ELISA Kit	78.13--5000 ng/mL	46.88 ng/mL	Sandwich
E-EL-H0470c	Human ApoE(Apolipoprotein E) ELISA Kit	23.44--1500 ng/mL	14.06 ng/mL	Sandwich
E-EL-H0488c	Human AQP-2(Aquaporin 2) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2500c	Human AR(Aldose reductase) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0043c	Human CRP(C-Reactive Protein) ELISA Kit	0.39--25 ng/mL	0.23 ng/mL	Sandwich
E-EL-H2069c	Human CTLA4(Cytotoxic T-Lymphocyte Associated Antigen 4) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0804c	Human C3(Complement 3) ELISA Kit	7.81--500 ng/mL	4.69 ng/mL	Sandwich
E-EL-H0058c	Human DPP4(Dipeptidyl Peptidase IV) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H3640c	Human EPO(Erythropoietin) ELISA Kit	1.56--100 mIU/mL	0.94 mIU/mL	Sandwich
E-EL-H1818c	Human GCK(Glucokinase) ELISA Kit	0.78--50 ng/mL	0.47 ng/mL	Sandwich
E-EL-H1917c	Human GHR/GHBP(Growth Hormone Binding Protein) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich

Cat. No.	Product Name	Detection Range	Sensitivity	Test Method
E-EL-H1919c	Human GHRL(Ghrelin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Competitive
E-EL-H2061c	Human GIP(Gastric Inhibitory Polypeptide) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H1491c	Human GJβ1(Gap Junction Protein Beta 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0148c	Human GLP-1(Glucagon Like Peptide 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H1825c	Human GLUT4(Glucose Transporter 4) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H2172c	Human HO1(Heme Oxygenase 1) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Sandwich
E-EL-H0086c	Human IGF-1(Insulin-like Growth Factor 1) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H0442	Human IGFBP-1(Insulin-like Growth Factor Binding Protein 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0103c	Human IL-10(Interleukin 10) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0149c	Human IL-1β(Interleukin 1 Beta) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H1030c	Human IL-2Rα(Interleukin 2 Receptor Alpha) ELISA Kit	93.75--6000 pg/mL	56.25 pg/mL	Sandwich
E-EL-H0102c	Human IL-6(Interleukin 6) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H0253c	Human IL-18(Interleukin 18) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H2665c	Human INS(Insulin) ELISA Kit	0.78--50 μIU/mL	0.47 μIU/mL	Sandwich
E-EL-H0452c	Human INSR(Insulin Receptor) ELISA Kit	1.56--100 ng/mL	0.94 ng/mL	Sandwich
E-EL-H5554c	Human IRS1(Insulin Receptor Substrate 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0113c	Human LEP(Leptin) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0115c	Human MAU(Microalbuminuria) ELISA Kit	3.13--200 μg/mL	1.88 μg/mL	Sandwich
E-EL-H0020c	Human MCP-1(Monocyte Chemotactic Protein 1) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0844c	Human MTHFR(Methylenetetrahydrofolate Reductase) ELISA Kit	78.13--5000 pg/mL	46.88 pg/mL	Sandwich
E-EL-H0096c	Human NGAL(Neutrophil Gelatinase Associated Lipocalin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1205c	Human NGF(Nerve Growth Factor) ELISA Kit	15.63--1000 pg/mL	9.38 pg/mL	Sandwich
E-EL-H0755c	Human NOS3/eNOS(Nitric Oxide Synthase 3, Endothelial) ELISA Kit	62.5--4000 pg/mL	37.5 pg/mL	Sandwich
E-EL-H1898c	Human NTRK2(Neurotrophic Tyrosine Kinase Receptor Type 2) ELISA Kit	93.75--6000 pg/mL	56.25 pg/mL	Sandwich
E-EL-H1014c	Human POMC(Pro-Opiomelanocortin) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H2298c	Human PON1(Paraoxonase 1) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1361c	Human PPAR-γ(Peroxisome Proliferator Activated Receptor Gamma) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1359c	Human PPARγC1α(Peroxisome Proliferator Activated Receptor Gamma Coactivator 1 Alpha) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H1237c	Human PYY(Peptide YY) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H1213c	Human RETN(Resistin) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H2183c	Human SAA(Serum Amyloid A) ELISA Kit	1.25--80 ng/mL	0.75 ng/mL	Sandwich
E-EL-H2115c	Human SOD2(Superoxide Dismutase 2, Mitochondrial) ELISA Kit	0.31--20 ng/mL	0.19 ng/mL	Competitive
E-EL-H0110c	Human TGF-β1(Transforming Growth Factor Beta 1) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich
E-EL-H0109c	Human TNF-α(Tumor Necrosis Factor Alpha) ELISA Kit	7.81--500 pg/mL	4.69 pg/mL	Sandwich
E-EL-H1511c	Human UCP-2(Uncoupling Protein 2, Mitochondrial) ELISA Kit	0.16--10 ng/mL	0.09 ng/mL	Sandwich
E-EL-H0111c	Human VEGF-A(Vascular Endothelial Cell Growth Factor A) ELISA Kit	31.25--2000 pg/mL	18.75 pg/mL	Sandwich



网 址：www.elabscience.cn

电 话：027-87526315

E-mail: techsupport@elabscience.cn tech@elabscience.cn

地 址：武汉市东湖新技术开发区高新大道 858 号生物医药园二期 B18 栋

